

BionicLoop Instructions for Use (IFU)

Document ID	IFU-BL-001
Version	1.17 (Approved for IDE submission)
Date	2026-09-04
Prepared by	Software Developer
Reviewed by	Camille Powe, MD, and Tara Bresnahan, RN; controlled correspondence dated 2026-09-03 and 2026-09-04

Investigational Use Only

Clinical Use Limitation: BionicLoop is investigational software and is not for independent clinical treatment use.

This document is intended for study-team review and supervised use preparation. The app screens shown are representative simulator captures from BionicLoop 1.0 (843). The alert and meal-availability reference tables in Part 2 are generated directly from the same controlled source and match the installed software.

Revision History

Version	Date	Author	Summary of Changes
1.0	2026-03-05	Team	Initial IFU draft package
1.1	2026-03-05	Team	Expanded operational narrative and alerts
1.2	2026-03-05	Team	Added simulator screenshots and workflow cleanup
1.3	2026-03-05	BionicLoop engineering	Added controlled title page, revision history, generated TOC, and print packaging
1.4	2026-04-05	BionicLoop engineering	Rebuilt as a workflow-driven user guide with refreshed screen set, operational checklists, alert-response guidance, and stronger task coverage
1.5	2026-07-09	BionicLoop engineering	Split into a participant Daily Use Guide and a Clinician & Study-Team Reference; refreshed every capture from the current software; rewrote Manual BG (digit keypad with confirmation tap), meal announcement (guided first step, carb comparison row, in-place timing recovery), and cancel-delivery flows; added Home alert-stack presentation, expanded Recent Dose Steps with status meanings, chart reading guide, and app-generated alert / meal-availability reference tables; added review sign-off page
1.6	2026-07-10	BionicLoop engineering	Added first-launch coverage: new section 10.3 (one-time Device Setup on a new phone, completed before sign-in, with a caution that the target profile starts on Pregnancy), new section 11.2 Clinical Unlock Setup (automatic and setup-code paths, one-time-per-phone behavior, single-use codes, retry pause), sign-in section now distinguishes a set-up phone from a new phone, target-change approval step documented, four new figures; figure index renumbered
1.7	2026-07-14	BionicLoop engineering	Added section 4.1 (If Sensor Readings Stop): fingerstick cadence during a sensor outage (every 2 hours; hourly after a reading over 200 or under 80), the Home countdown line, what Backup Basal Running means, and the due-soon/overdue alert tiers; alert reference tables regenerate from the app source and include the three new outage alerts
1.8	2026-07-14	BionicLoop engineering	Clinical review feedback (C. Powe): first-launch defaults are now Standard profile, 120 mg/dL target, 40-minute TMAX (Fiasp); the first-launch caution flips accordingly (select Pregnancy for participants assigned to the Pregnancy configuration); TMAX described as insulin action time; new section 11.5 Switching to a New Participant (New Participant Reset erases all participant data and sign-in for reused phones; the same-participant algorithm reset is relabeled); app install/reinstall note added; wizard and Clinical Settings figures recaptured; first-launch step regrouped into labeled Participant / Mode & Glucose Target / Insulin Settings boxes with per-mode target defaults (Standard 120, Pregnancy 110) and the review box now shows the mode; New Participant Reset offers to stop a running session first; empty-store loads return first-launch defaults without persisting them; Home chart figures recaptured with the demo live reading aligned to the chart series (no vertical jump at Now); notification-permission guidance added (readiness-checklist item and Section 7 note, 2026-07-14 audit)

Version	Date	Author	Summary of Changes
1.9	2026-07-24	BionicLoop engineering	Replaced the in-place second-tap Manual BG confirmation with a separate exact-value review alert; Change returns to the entered value and Use confirms it. Direct fingerstick-request states now share the Home gold BG-button emphasis. Added Part 1 Section 7, Using a Temporary Target (Pregnancy Mode): start/review/confirm flow with no preselection, active-state display on Home and the Exercise and Insulin Delivery screen, early end, exact-time automatic end with the Temporary Target Ended notification, activation-block messages, and a study-team note on primary-vs-safety/backup target separation (RA-022, clinical decision 2026-07-24); prior Sections 7-8 renumbered to 8-9. Settings walkthrough lists the new Exercise and Insulin Delivery screen (Temporary Target, Pod Activity suspension guidance, reset-location note). Clinical Settings documents the Standard-to-Pregnancy 90% meal-split draft default (stored values unchanged until Save). Sensor-outage section documents the overdue-fingerstick repeat (about every 15 minutes) and the countdown line flipping to overdue. Recent Dose Steps documents fallback-schedule rows (armed / schedule updated / schedule checked - unchanged, with per-segment rate provenance) and canceled/partial evidence persistence across relaunch per the RA-017 controlled clarification. Alert and meal-availability reference tables regenerated from app source. Three Temporary Target figures pending capture; figure numbers assigned at capture.
1.10	2026-07-30	BionicLoop engineering	Moved the existing Pod suspend/resume control directly onto Exercise and Insulin Delivery and removed its duplicate from normal Pod Settings. Clarified that the selected duration is a reminder only, insulin remains suspended until manually resumed, and suspend alerts open the relocated control.
1.11	2026-07-30	BionicLoop engineering	Moved Exercise and Insulin Delivery above User in Settings, aligned Temporary Target selectors/actions with the app control style, removed the unrelated participant-reset pointer, and added an explicit warning before target-only or duration-only selections are discarded.
1.12	2026-07-31	BionicLoop engineering	Added Dexcom G7 replacement-sensor recovery instructions for the always-visible Scan for new sensor action, persistent replacement progress, and the ten-minute Sensor Not Connected guidance. Clarified that reading age alone never changes the bound sensor.
1.13	2026-08-20	BionicLoop engineering	Clarified that an outlined red drop on the Home chart means a submitted fingerstick BG is waiting for a dosing step, and that the drop fills when the step uses it.
1.14	2026-08-21	BionicLoop engineering	Added deterministic Build 843 captures for Temporary Target selection, confirmation, and the active Home state; refreshed affected Build 843 workflow captures and completed the figure index.
1.15	2026-08-27	BionicLoop engineering	Aligned replacement-sensor instructions with RA-023: keep the study phone next to the intended replacement sensor during acquisition and verify the connected sensor/readings against the intended Dexcom app session before relying on BionicLoop.
1.16	2026-09-01	BionicLoop engineering	Recaptured the complete Build 843 figure set under controlled simulator conditions so all Home-screen figures show the intended FOR INVESTIGATIONAL USE ONLY banner and the exact v1.0 (843) build identity. No product behavior or shipping source changed.
1.16 approval	2026-09-03	Camille Powe, MD	Approved IFU-BL-001 revision 1.16 for IDE submission by controlled correspondence.
1.17	2026-09-04	Software Developer	Applied the approved labeling corrections: clarified Pregnancy configuration terminology, corrected Clinical Unlock section references, replaced an internal build-tool instruction, used plain punctuation, and documented the built-in investigational badge shown when no remote banner configuration is available. No product behavior or shipping source changed.

Table of Contents

- [Revision History](#)
- [Table of Contents](#)
- [Investigational Use Statement](#)
- [How This Guide Is Organized](#)
- [Part 1 - Daily Use Guide](#)
 - [1. The Home Screen](#)
 - [1.1 What Each Area Shows](#)
 - [1.2 When Alerts Are Showing](#)
 - [1.3 Before You Act](#)
 - [2. Reading the Chart](#)
 - [3. Signing In](#)
 - [4. Entering a Fingerstick BG](#)
 - [4.1 If Sensor Readings Stop](#)
 - [5. Announcing a Meal](#)
 - [6. Stopping a Meal Delivery](#)
 - [7. Using a Temporary Target \(Pregnancy Mode\)](#)
 - [7.1 Start a Temporary Target](#)
 - [7.2 While a Temporary Target Is Active](#)
 - [7.3 How a Temporary Target Ends](#)
 - [7.4 Also on the Exercise and Insulin Delivery Screen](#)
 - [8. Responding to Alerts](#)
 - [9. When to Contact the Study Team](#)
- [Part 2 - Clinician & Study-Team Reference](#)
 - [10. System Overview](#)
 - [10.1 System Components](#)
 - [10.2 What BionicLoop Does](#)
 - [10.3 What BionicLoop Does Not Do](#)
 - [11. Before First Use](#)
 - [11.1 Training and Readiness Checklist](#)
 - [11.2 Daily Start-of-Shift Checks](#)
 - [11.3 First Launch on a New Phone \(One-Time Setup\)](#)
 - [12. Settings and Clinical Configuration](#)
 - [12.1 Open Settings](#)
 - [12.2 Clinical Unlock Setup \(One-Time per Phone\)](#)
 - [12.3 Change Clinical Settings](#)

- [12.4 Start Algorithm Session](#)
- [12.5 Switching to a New Participant](#)
- [13. Device Workflows](#)
 - [13.1 Dexcom G7 First-Time Setup](#)
 - [13.2 Review CGM Status](#)
 - [13.3 Pod Setup](#)
- [14. Reviewing Recent Dose Steps](#)
 - [14.1 What the Step Statuses Mean](#)
- [15. Alert Reference](#)
 - [15.1 Alert-Response Principles](#)
- [16. Meal Availability Reference](#)
- [17. Troubleshooting and Recovery](#)
 - [17.1 If the Screen State and Device State Do Not Match](#)
- [18. Documentation and Escalation](#)
- [Appendix A. Figure Index](#)
- [Appendix B. Abbreviations and Acronyms](#)
- [Appendix C. Daily Use Quick Checklist](#)
- [Review and Approval](#)

Table of contents is generated automatically from document headings in the HTML and PDF outputs.

Investigational Use Statement

BionicLoop is investigational software and is **not approved for independent clinical treatment use**.

- Use only under approved study protocol and supervision.
- Use only after protocol-specific training on Dexcom G7, Omnipod DASH, and BionicLoop workflows.
- Do not rely on BionicLoop as the sole basis for treatment decisions outside approved investigational procedures.
- Follow site escalation procedures for therapy-impacting alarms, workflow failures, device faults, and documentation.

How This Guide Is Organized

Part	Audience	Contents
Part 1 - Daily Use Guide	Participants and trained supervised operators	The Home screen, entering a fingerstick BG, keeping dosing running when sensor readings stop, announcing meals, stopping a meal delivery, responding to alerts, when to contact the study team
Part 2 - Clinician & Study-Team Reference	Clinicians and study staff	System overview, readiness checklists, first launch on a new phone, clinical settings and algorithm session control, device workflows, dose-history review, complete alert and meal-availability reference tables, troubleshooting, documentation and escalation

This guide does not replace:

- Dexcom G7 app instructions,
- Omnipod DASH setup and safety instructions,
- the study protocol,
- site-specific escalation procedures.

Part 1 - Daily Use Guide

1. The Home Screen

Home is the first screen you see and the place to check before doing anything else.




Figure 1. Home with the current glucose value, Pod status, glucose and insulin chart, and the two action buttons.

1.1 What Each Area Shows

Home area	What it tells you
Glucose card (top left)	The most recent CGM value and trend arrow
Pod card	Insulin remaining in the Pod
Status ring (top right)	Green ring: automated dosing is running normally
Investigational-use badge	Shows the study-configured label. If no banner configuration is available, Home shows the built-in NOT FOR HUMAN USE label.
Alert area (above the chart)	The most important active alert, if any
Chart	Recent glucose (dots) and insulin given (bars), with 4h–24h views
Droplet button	Enter a fingerstick BG
Let's Eat button	Announce a meal

1.2 When Alerts Are Showing

When something needs attention, the most important alert appears as a card above the chart. If more alerts are active, a **"2 More"** control appears under the card - tap it to see the rest, and **"Show Less"** to collapse them again. The bell at the top shows the total count and opens the full Alert Center. When a related alert is hidden because the visible alert already explains it, a **"1 related alert in Alert Center"** link appears under the card - tap it to see everything.




Figure 2. Home showing the most important alert with more alerts collapsed behind it.




Figure 3. The same alerts expanded in place.

1.3 Before You Act

Do not announce a meal or enter a BG without checking first that:

- a current glucose value is showing (not missing or old),
- the Pod card shows an active Pod,
- no red (safety-critical) alert is unresolved.

If Home does not look ready, read the alert guidance in section 8 or contact the study team.

2. Reading the Chart

The chart shows glucose dots on top and insulin bars along the bottom. A few visual details carry meaning:

What you see	What it means
Solid colored bar	Insulin was given at that step
Bar with a yellow outline extending above it	Less insulin was given than planned - the outline shows the planned amount
Lighter, faded bars	Catch-up records added after a connection gap; the insulin was already accounted for
Small tick at the baseline	The app ran on schedule and no insulin was needed
Nothing at all	No dosing step ran at that time
Outlined red drop	A fingerstick BG was submitted and is waiting for a dosing step
Filled red drop	A dosing step used the submitted fingerstick BG

Touch and hold the chart to scrub: readouts appear above the chart showing the glucose and insulin values at the time under your finger, and disappear when you let go.

3. Signing In

On a phone that is already set up for the study, launching the app goes straight to sign-in. (The very first launch on a new phone instead asks for clinical settings before sign-in - that one-time setup is a study-team task, described in Part 2, Section 11.3.)

1. Launch the app.
2. Enter the account identifier provided by the study team.
3. Tap **Secure Sign In**.
4. Use **Create account** or **Forgot password?** only when the study team directs you to.
5. If you cannot sign in, contact the study team - do not keep retrying.




Figure 4. Sign-in screen.

If the app asks you to log in again while you are wearing the Pod, insulin delivery is not interrupted – signing back in restores remote monitoring for the study team. Sign back in when you can.

4. Entering a Fingerstick BG

Use the droplet button when the study workflow calls for a fingerstick blood glucose value. The value may be used for dosing decisions. It does **not** calibrate the CGM sensor.

1. On Home, tap the red droplet button.
2. Type the value on the keypad. The entry starts empty - the full value must be typed.
3. Tap Review BG.
4. Check the exact value in the review alert (for example, Use 143 mg/dL?). Tap Use 143 mg/dL to submit, or tap Change to return to the keypad without clearing the value.
5. Tap Cancel at any time to close without submitting.




Figure 5. Manual BG entry - the value starts empty.




Figure 6. A typed value ready to use.




Figure 7. The separate review alert shows the exact value before it is submitted.

Accepted values are 20–600 mg/dL. A value outside that range shows a red message and cannot be submitted.




Figure 8. An out-of-range value cannot be submitted.

4.1 If Sensor Readings Stop

If the sensor stops sending readings, the app keeps automated dosing running as long as you enter fingerstick BGs on schedule:

1. Home shows a line telling you when the next fingerstick is due ("Next fingerstick due by ..."). Enter a BG by that time using the red droplet button. Once the due time passes, the same line shows the check as overdue.
2. Checks are due every 2 hours. After a reading over 200 or under 80 mg/dL, the next check is due in 1 hour instead.

3. You will get a reminder about 15 minutes before a check is due. If a check becomes overdue, you will get an urgent alert, and it repeats about every 15 minutes until a BG is entered - enter one as soon as you can, or contact the study team.
4. If more than about an hour passes without any reading, the Pod switches to its steady backup basal rate and Home shows "Backup Basal Running." This is expected - entering a fingerstick BG puts automated dosing back in charge right away.
5. When sensor readings return, everything clears on its own and no action is needed.

There is no time limit: automated dosing continues as long as scheduled BGs keep coming in.

5. Announcing a Meal

Announce a meal just before eating, when Home is ready.

1. On Home, tap **Let 's Eat**.
2. Choose the meal: Breakfast, Lunch, or Dinner. The meal buttons pulse gently until you choose - that is the first required step.
3. Choose how the meal's carbs compare to your usual for that meal: **Less**, **Typical**, or **More**. **Typical** is preselected.
4. When the announcement is correct, slide the delivery control to deliver.




Figure 9. Meal announcement before a meal is chosen.




Figure 10. After choosing the meal, the slide-to-deliver control is enabled.

While the sheet is open you may briefly see **"Timing updated - a loop step just ran."** This is normal: the app re-checked its timing and kept your selections. Continue where you left off.

If meal announcement is not available right now, the app says why (for example, the Pod is still delivering, or a recent delivery has not been confirmed yet) and, when the reason is time-based, when to try again. The full list of these messages is in Part 2, section 16.

6. Stopping a Meal Delivery

If a meal announcement is delivering and needs to stop, the meal sheet shows the delivery in progress with a red **Cancel Delivery** slider.

1. Tap **Let's Eat** if the meal sheet is not already open. While a meal delivery is in progress it opens to the delivery-in-progress view.
2. Confirm this is the delivery you intend to stop.
3. Slide **Cancel Delivery**.
4. Wait for Home to show the summary of what was actually delivered.




Figure 11. A meal delivery in progress with the cancel slider.

After cancelling, Home shows an orange summary above the chart with the amount requested and the amount actually delivered before the stop.




Figure 12. The partial-delivery summary after a cancelled meal delivery.

What the summary means:

- the delivered amount is what reached the Pod before the stop,
- that insulin stays part of the recorded history and is accounted for in later dosing,
- the summary stays visible long enough to review and then clears on its own.

If a meal announcement is still needed after a cancellation, open **Let's Eat** again and enter it deliberately.




Figure 13. Announcing a new meal after a cancelled delivery.

7. Using a Temporary Target (Pregnancy Mode)

A Temporary Target raises your glucose target for a short, fixed time - for example around exercise - without changing your regular settings. It is available only when the phone is configured in Pregnancy mode; in Standard mode the screen shows "Temporary Target is available in Pregnancy mode." and no selection is offered.

7.1 Start a Temporary Target

1. Open Settings (the gear on Home), then **Exercise** and **Insulin Delivery**.

2. Under **Temporary Target**, choose a target: 120 mg/dL or 130 mg/dL.
3. Choose a duration: 30 min, 1 hr, 1 hr 30 min, or 2 hr. Nothing is preselected - both choices are always yours to make.
4. Tap **Review Temporary Target**. A confirmation sheet shows the target, the duration, and the exact end time.
5. Tap **Start** to begin, or **Cancel** to go back without changing anything.

If you select either a target or duration and then try to leave this screen before starting, the app warns that the selections will be lost. Tap **Keep Editing** to return, or **Discard Selections** to leave without starting a Temporary Target.

The screenshot shows a mobile app interface for 'Exercise and Insulin Delivery'. At the top, the status bar displays the time 9:41, signal strength, Wi-Fi, and battery icons. The app title 'Exercise and Insulin Delivery' is centered at the top of the screen. Below the title, the section 'Temporary Target' is displayed. Under this section, there are two sub-sections: 'Target' and 'Duration'. The 'Target' section has two buttons: '120 mg/dL' and '130 mg/dL ✓', with the latter being selected. The 'Duration' section has three buttons: '30 min', '1 hr ✓', and '2 hr', with '1 hr' being selected. Below these sections is a blue button labeled 'Review Temporary Target'. At the bottom of the screen, there is a section titled 'Suspend Insulin Delivery' which includes a pause icon and the text 'Suspend Insulin Delivery'. Below this, a message states: 'An active Pod is required to suspend or resume insulin delivery.'

Figure 14. Temporary Target selection: target and duration.

9:41

< Exercise and Insulin Delivery

Temporary Target

Target

120 mg/dL 130 mg/dL ✓

Duration

30 min 1 hr ✓

Confirm Temporary Target

Target 130 mg/dL

Duration 1 hr

Ends 9:38 AM

Cancel Start

Figure 15. Confirmation sheet with target, duration, and end time.

If the app cannot start a Temporary Target, it says why:

Message	What it means
Start the algorithm before using Temporary Target.	The algorithm session is not running. Starting it is a study-team task (Part 2, Section 12.4).
Set up a Pod before using Temporary Target.	No active Pod (Part 2, Section 13.3).

Message	What it means
Wait for backup insulin recovery to finish, then try again.	The app is still reconciling backup-basal delivery after a connection gap. Wait for Home to settle, then retry.
Choose a target and duration, then try again.	One of the two selections is missing.

7.2 While a Temporary Target Is Active

- Home shows a status line with the target and end time - for example, "Temporary Target 120 mg/dL · until 3:45 PM".
- The **Exercise** and **Insulin Delivery** screen shows the **Target**, when it **Ends**, and "Temporary Target is active", with an **End Temporary Target** button.
- Automated dosing continues and aims at the temporary target. Safety checks and the Pod's backup basal schedule stay on your permanent settings (see the study-team note below).




Figure 16. Home status line while a Temporary Target is active.

7.3 How a Temporary Target Ends

- **On its own:** it ends exactly at the end time, even if the phone was locked or the app was closed in between. A notification says "Temporary Target Ended - Your usual glucose target is active again."
- **Early:** tap End Temporary Target on the same screen, then confirm ("Your usual glucose target will be active again").
- It also ends automatically if the algorithm session is stopped, the target profile is changed, or a participant reset is performed.

7.4 Also on the Exercise and Insulin Delivery Screen

- **Suspend Insulin Delivery:** choose a reminder for 30 minutes, 1 hour, 1 hour 30 minutes, or 2 hours. Insulin delivery stops and remains suspended until you return to this screen and tap **Resume Insulin Delivery**; the reminder does not resume insulin automatically.
- If no usable Pod is active, or the Pod is still changing its delivery state, the suspend/resume action is unavailable until the Pod is ready.

Note for the study team. Only the primary dosing target follows the Temporary Target, and it is resolved by each algorithm step's execution time: the first step at or after the end time uses the permanent target even if the end notification was never delivered. The safety controller and the Pod's programmed backup-basal schedule remain on the permanent Pregnancy configuration throughout (clinical decision 2026-07-24; risk analysis RA-022).

8. Responding to Alerts

1. Read the alert card on Home. Each alert says what is happening and what to do.
2. Handle the top alert first - it is the most important one and often explains the others.
3. Tap the bell to open Alert Center for the complete list, including recently cleared alerts.
4. Do not assume things are back to normal until Home shows the alert has cleared.




Figure 17. Alert Center with an active alert.

Glucose alarms come from the Dexcom app. BionicLoop's CGM messages are extra context inside the app only - BionicLoop does not sound alarms for low or high glucose. Keep the Dexcom G7 app's alarms configured as the study team directed.

Keep BionicLoop notifications on. Fingerstick-due and overdue reminders arrive as iPhone notifications when the app is in the background. If notifications are turned off for BionicLoop in iPhone Settings, those reminders will not appear. The study team confirms this permission during setup - tell them if it has been changed.

9. When to Contact the Study Team

Contact the study team rather than retrying when:

- a red (safety-critical) alert does not clear after following its guidance,
- the app and the physical devices disagree (for example, the app shows an active Pod but no Pod is worn),
- a delivery outcome is reported as unconfirmed and does not resolve,
- you cannot sign in and the app is needed for the study workflow,
- anything about dosing looks wrong to you - when in doubt, call.

Before calling, if possible: note the time, what you were doing, and take a screenshot of Home and the alert.

Part 2 - Clinician & Study-Team Reference

10. System Overview

BionicLoop is a research closed-loop app that uses Dexcom G7 glucose data, Omnipod DASH pump connectivity, and a 5-minute algorithm step cadence to run algorithm-directed insulin dosing under study supervision.

10.1 System Components

Component	Role in workflow
iPhone running BionicLoop	Home view, settings, meal/BG entry, alert presentation, telemetry
Dexcom G7 sensor and Dexcom ecosystem	Source of glucose data and primary CGM alarm handling
Omnipod DASH Pod	Insulin delivery path
BionicLoop cloud/telemetry path	Remote monitoring and quality review support when configured

10.2 What BionicLoop Does

- runs the dosing algorithm on a 5-minute step cadence during an active session,
- displays current CGM, Pod, and loop status on Home,
- supports clinical configuration and algorithm session control,
- accepts manual BG entry and meal announcements,
- records recent algorithm/dosing steps for review,
- presents in-app alerts and recent alert history.

10.3 What BionicLoop Does Not Do

- BionicLoop does **not** replace Dexcom G7 app alarming.
- BionicLoop does **not** send OS notifications for CGM alerts.
- BionicLoop CGM alerts are informational on Home / Alert Center.
- The Dexcom G7 app remains the **source of truth for CGM alarms** and should be configured per protocol.

CGM alerting policy: Keep the FDA-cleared Dexcom application configured correctly for urgent low, low, high, signal, and other sensor alarms. Review any BionicLoop CGM messages as supplemental app context, not as the primary alarm source.

11. Before First Use

11.1 Training and Readiness Checklist

1. The user has been trained on the study workflow and escalation rules.
2. The iPhone is charged and Bluetooth is enabled.
3. BionicLoop notifications are allowed on the iPhone (Settings > Notifications > BionicLoop).
Safety reminders - like a fingerstick coming due during a sensor outage - arrive as phone notifications and stay silent if this permission is off.
4. The user can sign in or site login support is available.
5. Dexcom G7 app access is available for pairing and sensor management.
6. Pod supplies and site procedures are available.
7. Subject ID, weight, and approved clinical settings are available. On a new phone, the one-time first-launch setup (Section 11.3) is complete or scheduled.
8. The user understands that unresolved therapy-impacting issues require escalation, not repeated blind retries.

11.2 Daily Start-of-Shift Checks


1. Launch BionicLoop and confirm the correct subject context.
2. Confirm CGM value freshness and trend availability.
3. Confirm Pod presence and connection state.
4. Review active alerts before announcing a meal or entering BG.
5. If using remote monitoring, confirm login/session continuity.

11.3 First Launch on a New Phone (One-Time Setup)

The very first time BionicLoop opens on a new phone, it shows a two-step `Device Setup` screen **before sign-in**. This is a study-team task, done once per phone. The app cannot reach the Home screen until this setup and the first sign-in are complete.

1. Launch BionicLoop. The `Device Setup` screen appears (Step 1 of 2 - Enter clinical settings).
2. In `Participant`, enter the Subject ID and weight (in pounds) from the approved subject plan.
3. In `Mode & Glucose Target`, pick the mode: **starts on Standard (target defaults to 120 mg/dL)**. Selecting `Pregnancy` switches the target options and defaults to 110. Change the target only if the approved subject plan says so.

4. In **Insulin Settings**, set the meal insulin given up front (75% or 90%) and TMAX - the insulin action time (defaults to 40 minutes for Fiasp).
5. Check the **Review** box - it repeats the mode and every value you are about to save.
6. Tap **Save Settings**.



9:41

Device Setup

Complete initial clinical settings before signing in to BionicScout.

● ● Step 1 of 2

Enter clinical settings

Study staff enter these values from the study protocol before sign-in.

Participant

Subject ID

Subject ID

Weight (pounds)

Weight in pounds

Mode & Glucose Target

Mode

Pregnancy Standard

Glucose target

110 120 130

mg/dL. Standard defaults to 120; Pregnancy defaults to 110.

Insulin Settings

Meal insulin given up front

75% 90%

TMAX - insulin action time

40 min ↕

Figure 18. First launch: enter clinical settings before sign-in.

7. The app confirms the setup is saved and moves to **Step 2 of 2**. Tap **Sign in** to BionicScout and sign in (the sign-in screen is shown in Part 1, Section 3).




Figure 19. First launch: setup saved, ready to sign in.

8. After sign-in, finish first-time setup from the Home settings gear:

- Dexcom G7 setup (Section 13.1).
- Pod setup (Section 13.3).
- Clinical Unlock Setup (Section 12.2), so settings can be changed later.

Notes:

- This first save does not ask for an unlock code. Every change after it requires Clinical Unlock (Sections 12.2 and 12.3).

- Double-check the target profile before saving. If a participant assigned to the Pregnancy configuration is left on the Standard default, the algorithm will use higher glucose targets than intended for pregnancy. Fixing it afterward requires Clinical Unlock.
- Installing the app: the study team provides the install link. Reinstalling after a deletion uses the same link; a reinstalled phone starts at this first-launch setup.

12. Settings and Clinical Configuration

The screenshot shows the 'Bionic Loop Settings' screen of a mobile application. At the top, the status bar displays the time 9:41, signal strength, Wi-Fi, and battery icons. The app title 'Bionic Loop Settings' is centered, with a 'Done' button in the top right corner. Below the title is a 'Clinical Settings' section with a right-pointing arrow. This section contains several settings: 'Target Access Profile' set to 'Pregnancy', 'Glucose Target [mg/dL]' with three buttons (90, 100, 110) where 110 is selected, a note 'Changing target requires clinical approval capture.', 'Subject ID' set to 'UITEST', 'Weight (lb)' set to 150, 'Meal Split [% upfront]' set to 75%, and 'TMAX [minutes]' set to 40. Below this is an 'Exercise and Insulin Delivery' section with a running person icon and a right-pointing arrow. The 'User' section at the bottom includes 'Account & Session' set to 'Signed In' and 'Recent Dose Steps' set to 11, both with right-pointing arrows.

Setting	Value
Target Access Profile	Pregnancy
Glucose Target [mg/dL]	110
Subject ID	UITEST
Weight (lb)	150
Meal Split [% upfront]	75%
TMAX [minutes]	40
Exercise and Insulin Delivery	
Account & Session	Signed In
Recent Dose Steps	11

Figure 20. Bionic Loop Settings sheet.

12.1 Open Settings

1. Tap the settings gear on Home.
2. Review the current target, subject ID, weight, and algorithm-input summary.
3. Open **Exercise and Insulin Delivery** for the Temporary Target (Pregnancy mode, Part 1 Section 7) and for insulin-suspension guidance.
4. Open device settings from this sheet when CGM or Pod setup/review is needed.

12.2 Clinical Unlock Setup (One-Time per Phone)

Changing saved clinical settings requires an unlock code, and unlock codes only work after a one-time setup that ties the phone to the subject. Do this once during first-time setup, after signing in.

1. Open Settings, then **Clinical Settings**. If the screen says "This device is not set up to accept unlock codes yet," complete this setup first.
2. Preferred path: confirm the Subject ID shown, then tap **Set Up Automatically**. This uses the app's BionicScout sign-in - the phone must be signed in, and the subject must be assigned to this phone (that happens automatically after the first sign-in).
3. Alternative path: tap **Have a setup code? Enter it manually**, enter the setup code provided by a study administrator from the BionicScout dashboard, then tap **Save Setup Code**.
4. When setup is complete, the screen shows the unlock setup as ready.




Figure 21. Clinical Unlock Setup before the one-time setup is completed.

About unlock codes:

- Each unlock code works once. Get a fresh code from the study team (BionicScout dashboard) each time you need to unlock.
- An unlock ends on its own after a short time, or when you tap **Lock Clinical Settings**.
- After 5 wrong codes, the screen pauses new attempts for about a minute.

12.3 Change Clinical Settings

Clinical Settings are unlock-code protected and should only be changed by authorized users.

1. Open Settings, then **Clinical Settings**.
2. Enter the clinical unlock code and tap **Unlock Clinical Settings**.
3. Review or update: Subject ID, Weight (lb), glucose target, meal split, TMAX. Switching the target profile from **Standard** to **Pregnancy** drafts the meal split to 90% (the **Pregnancy** default); nothing stored changes until **Save**.
4. Changing the regular glucose target also asks for the name of the clinician who approved the change.
5. Save only after confirming the values match the approved subject plan.




Figure 22. Clinical Settings locked: enter a one-time unlock code.

9:41

< Clinical Settings

Clinical Settings

Clinician controls unlocked.

[Lock Clinical Settings](#)

Subject

UITEST

Weight (lb) 150

Algorithm Inputs

Regular Settings Target Range:

Pregnancy Standard

Glucose Target [mg/dL]:

90 100 110

Clinical Settings and Regular Settings are both constrained to the Pregnancy range: 90 / 100 / 110 mg/dL.

Meal Split [% upfront]:

Figure 23. Clinical Settings unlocked, with algorithm input controls.

12.4 Start Algorithm Session

Use **Start Algo** only from Clinical Settings, and only when no active algorithm session is already running.

Before starting:

1. Confirm subject ID and weight.
2. Confirm the approved target profile and input settings.
3. Confirm CGM and Pod workflows are ready enough to support the first anchored loop step.

12.5 Switching to a New Participant

Study phones are reused. When a participant finishes the study, use `New Participant Reset` in `Clinical Settings` to prepare the phone for the next person. (The `Same-Participant Reset` button next to it is different - it restarts the algorithm for the same person and keeps their settings.)

Before the reset, deactivate and remove the Pod (Pod settings) - the reset refuses to run while a Pod is active.

Then:

1. Open `Settings`, then `Clinical Settings`, and unlock.
2. Tap `New Participant Reset`. If an algorithm session is still running, the app first offers `Stop Session and Continue` - confirm it.
3. Read the confirmation and tap `Erase and Start New Participant`.
4. The phone erases the participant's settings, dosing history, alerts, sensor and Pod pairing, unlock setup, and sign-in, and returns to first-time setup (Section 11.3) ready for the next participant.

13. Device Workflows

13.1 Dexcom G7 First-Time Setup

1. Open `Settings`, then `Open CGM Settings`.
2. Follow the guided Dexcom G7 setup path.
3. Continue pairing/calibration/sensor-management tasks in the Dexcom app as required.




Figure 24. Dexcom G7 setup entry view.

13.2 Review CGM Status

Check sensor expiration timing, last reading value and timestamp, trend information, and connection recency.

Scan for new sensor remains available even while BionicLoop is searching for the sensor already in use. Use it only after starting a replacement G7 or when the study team directs you to replace the current sensor. An old reading or a temporary connection gap does not make BionicLoop release the current sensor.

After you choose **Scan for new sensor**, a **Scanning for new sensor...** row shows that replacement acquisition is active. The replacement search survives closing and reopening BionicLoop. If the new sensor has not connected after 10 minutes, BionicLoop shows **Sensor Not Connected**. Open the Dexcom app to finish starting the replacement sensor, then return to BionicLoop Settings and try again. Do not select replacement scanning merely to recover a temporary signal gap from the sensor already in use.

Keep the study phone next to the intended replacement sensor and away from other active G7 sensors while BionicLoop acquires it. After BionicLoop connects, compare the sensor status and glucose reading with the intended session in the Dexcom app before relying on BionicLoop for automated dosing. Do not select a sensor by a Bluetooth name or suffix.





Figure 25. Dexcom G7 operational/settings view.

13.3 Pod Setup

1. Open Settings, then Open Pod Settings.

2. Follow the Pod setup instructions for reminder configuration, fill, pairing, and placement.
3. Complete all physical device steps per Omnipod instructions.




Before proceeding with Pod-dependent actions, confirm: the Pod is active, signal is present, no pump fault/incompatible/setup-incomplete alert is active, and the Pod is not beyond its usable service window.

14. Reviewing Recent Dose Steps

Use **Recent Dose Steps** (in Settings) to review what the algorithm recently asked for and what the Pod confirmed.

1. Open Settings, then **Recent Dose Steps**.
2. Each row shows the step, the insulin amount, the CGM value used, and a status.
3. Tap a row to expand its detail, including timestamps and the meal context when one applies.



Recent Dose Steps			
Step 472	Requested	0.05 U	8:32:20 AM >
Step 471	Requested	0.05 U	8:27:20 AM >
Step 470		0.00 U	8:22:20 AM >
Step 469	Requested	0.05 U	8:17:20 AM >
Step 468	Requested	0.80 U	8:12:20 AM >
Step 467		0.00 U	8:07:20 AM >
Step 466	Requested	0.05 U	8:02:20 AM >
Step 465		0.00 U	7:57:20 AM >
Step 464	Requested	0.05 U	7:52:20 AM >
Step 463		0.00 U	7:47:20 AM >

Figure 27. Recent Dose Steps.




Figure 28. An expanded step row with detail.

14.1 What the Step Statuses Mean

The step list distinguishes what was **requested** from what the Pod **confirmed**:

Status	Meaning
Delivered	The Pod confirmed the insulin was delivered
Delivering	The Pod is actively delivering this dose
Requested	The dose command was sent; Pod confirmation has not been received yet

Status	Meaning
Assumed delivered	Confirmation was impossible (for example, the Pod was replaced); the amount is treated as delivered per clinical policy
Replayed	A catch-up record written after a connection gap so the history stays complete
No dose	The step ran and no insulin was needed

Alongside step rows, `Fallback` schedule rows appear whenever the Pod's programmed backup schedule is armed, updated, or checked and found unchanged. Each shows the four programmed segment rates and, per segment, whether the rate was observed or imputed.

Canceled deliveries keep their evidence across an app relaunch. A partial that was capped while delivery was still active may later show a higher settled amount once the same Pod confirms it - other outcomes do not change on their own.

Treat a persistent `Requested` or an unexpected `Assumed delivered` as worth reviewing: confirm current Pod state and the surrounding steps.

15. Alert Reference

The table below is generated directly from the app's alert catalog and matches the installed software. Severity meanings: **Safety-critical** alerts affect active therapy and require action; **Action needed** alerts require a user response; **Informational** alerts provide context only.

On Home, when several alerts are active at once, the most severe shows first and generic alerts are hidden while a more specific active alert already explains the same problem. When that happens, a tappable "related alert in Alert Center" line appears under the visible card, so the bell count and the Home surface never silently disagree. Alert Center always shows the complete list.

Alert	Severity	What it means	What to do
Log In Required	Action needed	Remote monitoring is unavailable. Local therapy remains available on this phone.	Log in now to restore remote monitoring.
Subject ID Conflict	Action needed	Subject ID already in use.	Choose a different subject ID and save again.
Clock Check Needed	Action needed	Your phone time appears incorrect.	Check Date & Time settings and enable Set Automatically.

Alert	Severity	What it means	What to do
Pump Signal Loss	Action needed	Pump communication is unavailable.	Keep phone near pod and verify Bluetooth.
No Active Pod	Safety-critical	Current pod is inactive or deactivated.	Start Pod setup and pair a new pod.
Low Reservoir	Action needed	Pod insulin volume is low.	Prepare a replacement pod soon.
Pod Expiring	Action needed	Pod expires in 30 minutes.	Replace the Pod within 30 minutes to avoid missed insulin.
Pod Expired	Safety-critical	Pod is expired.	Replace pod now before grace period ends.
Pump Fault	Safety-critical	Pod reported a fault and delivery may be interrupted.	Replace pod and follow recovery workflow.
Incompatible Pod	Safety-critical	Connected pod is not compatible with this system.	Pair a compatible pod and complete setup again.
Pod Setup Incomplete	Action needed	Pod setup is not complete.	Complete pod setup flow.
Delivery Suspended	Action needed	Insulin delivery is currently suspended.	Open Settings > Exercise and Insulin Delivery to resume insulin delivery when appropriate.

Alert	Severity	What it means	What to do
Resume Insulin	Action needed	The suspension reminder has ended. Insulin delivery remains suspended.	Tap Resume Insulin to restart delivery.
Insulin Still Suspended	Safety-critical	Insulin delivery has remained suspended for at least 15 minutes after the reminder ended.	Tap Resume Insulin now. If delivery does not resume, keep the phone near the Pod and try again.
Pump Time Mismatch	Action needed	Pump time differs from device time.	Review and sync pump time in Pump settings.
CGM Unavailable	Informational	Reliable CGM data is not currently available.	Check sensor status and wait for reliable readings.
Urgent Low Glucose	Safety-critical	Current glucose is below 55 mg/dL.	Confirm participant safety and treat low glucose per protocol.
Algorithm Stepping Interrupted	Action needed	Automated dosing has not run in the last 15 minutes.	Open the app and check the Home screen alerts.
Dosing Paused to Re-Sync	Informational	The app is re-syncing its step count with the algorithm.	No action needed. Dosing resumes automatically within a few minutes.
CGM Failed/Expired	Informational	Sensor is failed or expired.	Replace sensor and restore CGM coverage.
Awaiting First Step	Action needed	Waiting for the first glucose reading to begin automated dosing.	Wait for next valid CGM reading.
Loop Command Blocked	Action needed	Dosing was calculated but can't reach the Pod.	Move the phone closer to the Pod and check Bluetooth.

Alert	Severity	What it means	What to do
Meal Not Available	Informational	Meal announce is unavailable at this time.	Try again at the time shown.
Blood Sugar Check Requested	Action needed	A fingerstick BG check or sensor read is needed for automated dosing to resume.	Tap the red droplet button and enter the value.
Dosing Paused	Informational	The dosing calculator is waiting for glucose readings. Dosing resumes automatically when readings return, or sooner if you enter a fingerstick BG.	No action is needed. If this lasts more than 15 minutes, check your sensor.
Dosing Paused - Sensor Data Needed	Safety-critical	Automated dosing has been paused for over 15 minutes waiting for glucose readings. Dosing resumes automatically when sensor readings return, or sooner if you enter a fingerstick BG with the red droplet button.	Check your sensor and its connection. Contact the study team if readings do not return.
Backup Insulin Not Confirmed	Safety-critical	About 0.35 U of backup insulin given during a connection gap could not be confirmed. Contact the study team.	Review Recent Dose Steps and contact the study team if this was unexpected.
Fingerstick BG Due Soon	Action needed	Enter a fingerstick BG in the next 15 minutes to keep automated dosing running. Not needed if sensor readings return first.	Tap the red droplet button and enter the value.
Backup Basal Running	Informational	No glucose readings for over an hour. The pod is delivering its backup basal rate. Dosing resumes automatically when sensor readings return, or sooner if you enter a fingerstick BG.	Enter a fingerstick BG when you can, or wait for the sensor to reconnect.
Fingerstick BG Overdue	Safety-critical	The pod is delivering only its backup basal rate. Dosing resumes automatically when sensor readings return, or sooner if you enter a fingerstick BG with the red droplet button.	Enter a fingerstick BG now, or contact the study team.

15.1 Alert-Response Principles

1. Check whether the alert changes therapy readiness.

2. Confirm the device state directly when applicable.
3. Review Home and recent dose context.
4. Escalate unresolved therapy-impacting conditions instead of repeatedly retrying actions.

16. Meal Availability Reference

When meal announcement is unavailable, the composer explains why with one of the following messages (generated directly from the app source). Time-dependent messages show the actual time in the app; {time} marks where it appears.

Message shown	Offers Open Settings	Condition
Start Bionic Loop to meal announce.	Yes	loopOff
Enter Subject ID and Weight in Settings.	Yes	missingProfile
Connect a Pod before meal announce.	No	noPump
Backup insulin may be active after a connection gap. Reconnect the Pod and wait for recovery to finish.	No	fallbackReconciliationRequired
Meal announce is available after the first loop step.	No	awaitingFirstStep
Meal announce opens at {time}.	No	tooEarly
Meal announce is unavailable for the current step. Try again at {time}.	No	tooLate
Another loop step executed while meal announce was open. Review current status and re-enter the meal if still needed.	No	slotConflict
Pump is delivering. Try again when delivery completes.	No	pumpDelivering
Insulin delivery is paused. Resume insulin delivery before announcing a meal.	No	pumpSuspended
Waiting for a recent CGM value.	No	missingFreshGlucose
Signal Loss. Reconnect Pod and try again.	No	signalLoss
Meal announce already in progress.	No	requestInProgress

Message shown	Offers Open Settings	Condition
A previous meal announcement has an uncertain delivery outcome. Verify insulin delivery before announcing another meal.	No	uncertainPreviousMeal
A recent insulin delivery hasn't been confirmed yet. Reconnect the Pod before announcing a meal.	No	unreconciledIssuedDose
The app is re-syncing its step count. Meals are available again once dosing resumes in a few minutes.	No	algorithmStepDrift
Meal was not announced. Backup basal is running because glucose readings are unavailable. Enter a fingerstick BG or wait for a sensor reading, then try again. No meal insulin was delivered.	No	backupBasalRunning
Meal announcement is unavailable right now.	No	generic fallback

If the composer is already open when a condition arises, the message renders inside the composer with the user's selections preserved, and clears on its own if the condition resolves.

17. Troubleshooting and Recovery

Situation	What to do
No current CGM value shown	Check the Dexcom app, sensor state, and reading freshness before relying on BionicLoop glucose context
Meal announcement unavailable	Read the message and check section 16; confirm Pod, CGM freshness, timing window, and session status
Pod unavailable / no active Pod	Reconnect or replace the Pod per protocol before proceeding
Pump fault / incompatible alert	Treat as therapy-impacting; replace/recover per protocol and verify Home clears appropriately
Pod expired or stopping soon	Plan Pod replacement before continued therapy reliance
Partial delivery recorded	Review the delivered amount, confirm whether further meal announcement is still needed, and document if clinically relevant

Situation	What to do
Backup insulin during a connection gap could not be confirmed	Contact the study team; review Recent Dose Steps and the surrounding replayed records
Login required for remote monitoring	Reauthenticate when allowed; document the monitoring gap if required

17.1 If the Screen State and Device State Do Not Match

1. Stop and reassess.
2. Review active alerts.
3. Verify Dexcom state in the Dexcom app and Pod state in the Pod workflow.
4. Capture screenshots if needed.
5. Escalate per site process.

18. Documentation and Escalation

For unresolved or therapy-impacting conditions:

1. Stabilize the situation per study protocol.
2. Record the time and sequence of events.
3. Capture the relevant Home / Alert Center / device screen.
4. Note any partial delivery amount or relevant recent step context.
5. Escalate to the supervising clinical/study contact.

Appendix A. Figure Index

Figure	File	Section use
Figure 1	home-default.png	Home orientation
Figure 2	home-alert-stack.png	Alert stack, collapsed
Figure 3	home-alert-stack-expanded.png	Alert stack, expanded
Figure 4	auth-login-start.png	Sign in
Figure 5	home-manual-bg.png	Manual BG, empty entry
Figure 6	manual-bg-entered.png	Manual BG, value typed
Figure 7	manual-bg-confirm.png	Manual BG, exact-value review alert
Figure 8	manual-bg-out-of-range.png	Manual BG, out-of-range
Figure 9	home-meal-announcement.png	Meal composer
Figure 10	meal-ready-to-deliver.png	Meal ready-to-deliver
Figure 11	home-meal-cancel-delivery.png	Cancel-delivery control
Figure 12	home-meal-cancel-delivery-partial.png	Partial-delivery summary
Figure 13	home-meal-cancel-delivery-ready.png	Re-announce after cancel
Figure 14	temporary-target-selection.png	Temporary Target selection
Figure 15	temporary-target-confirm.png	Temporary Target confirmation
Figure 16	home-temporary-target-active.png	Temporary Target active (Home)
Figure 17	home-alert-center.png	Alert Center
Figure 18	clinical-onboarding-settings.png	First launch, clinical settings step
Figure 19	clinical-onboarding-signin.png	First launch, ready to sign in
Figure 20	home-settings.png	Settings
Figure 21	clinical-unlock-setup.png	Clinical unlock one-time setup
Figure 22	clinical-unlock-locked.png	Clinical settings locked, unlock entry
Figure 23	home-clinical-settings.png	Clinical settings unlocked
Figure 24	cgm-setup.png	CGM setup

Figure	File	Section use
Figure 25	cgm-settings.png	CGM settings
Figure 26	pump-setup.png	Pod setup
Figure 27	recent-dose-steps.png	Recent dose review
Figure 28	recent-dose-steps-expanded.png	Recent dose detail

Appendix B. Abbreviations and Acronyms

Abbreviation	Meaning
IFU	Instructions for Use
CGM	Continuous Glucose Monitor
BG	Blood Glucose
TMAX	Time to maximum insulin effect parameter
DASH	Omnipod DASH pump system

Appendix C. Daily Use Quick Checklist

1. Confirm the correct subject and active session context.
2. Confirm CGM freshness.
3. Confirm Pod readiness.
4. Review alerts before acting - top alert first.
5. If a Temporary Target is active, Home shows it with its end time - confirm it matches the plan.
6. Use the droplet button or Let 's Eat only when Home is ready.
7. Review any partial delivery or interruption summary before proceeding.
8. Document and escalate any therapy-impacting mismatch or unresolved alert.

Review and Approval

This version was approved for IDE submission by controlled correspondence on 2026-09-04. The table below may be completed if the sponsor's document-control process requires additional signatures.

Role	Name	Signature	Date
Software Developer			
Clinical investigator or principal investigator			
Sponsor representative			

Approval of this IFU version covers the workflows, screens, and generated reference tables as of the document date. Controlled screenshots and reference tables must be regenerated from the approved application baseline before any subsequent revision.